PRODUKTINFORMATION

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

Verfügbar in:
3shape Model Builder
exocad model creator
BISS Promadent
SHERAeasy-model

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


Kategorie Champions Implants +

System Champions - Revolution V2.14.1

Kit Champions - Revolution V2.14.1
Champions - Revolution V2.14.2

REF 31219 DIM Analog
REF31131 DIM Analog Multi-Unit

Abutment


Champions Titan

Kategorie Champions Implants +

System Champions - Revolution V2.14.2

Kit Champions - Revolution V2.14.1
Champions - Revolution V2.14.2

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Champions DIM Model

CONTENT

1. Options for construction
2. Options for scanning Champions scanbodies
3. Installation/Instructions for selecting exocad
4. Installation/Instructions for selecting 3shape
5. Instructions for Laboratory Analogs for 3D-printed models

1. What options for construction are available in the libraries for 3shape and exocad?

With scanbody Impression Post Peek REF 3122 + Gingiva-Shuttle REF 3100
or Scan-Abutment Champions (R)Evolution® REF 31230

Preform C2 REF 31071

Custom abutment titanium
with Sirona® base

PreFace_V2_11,5 REF 31072
(Abutment Preform ME-S)

Custom titanium abutment
with Medentika® base

PreFace_V2_16 REF 31073
(Abutment Preform ME-L)

Custom titanium Abutment
with Medentika® base

Klebebasis_GH1 REF 3104
(Abutment Adhesive GH1)

Adhesive-based constructions:
e.g., screw-retained single crown, custom abutment
(zircon), primary crown (zircon, NEM [non-
precious metal alloy]), no screw-retained bridges!

Klebebasis_GH1_L REF 31043
(Abutment Adhesive L GH1)

Screw-retained single crown in the
anterior and premolar regions

Klebebasis_GH2 REF 31044
(Abutment Adhesive GH2)

Screw-retained single crown in the
anterior and premolar regions

Klebebasis_GH3 REF 31045
(Abutment Adhesive GH3)

Screw-retained single crown in the
anterior and premolar regions

ICA-Klebebasis_GH1 REF 3102
(ICA-Abutment GH1)

Screw-retained single crown in the molar region

ICA-Klebebasis_GH2 REF 31021
(ICA-Abutment GH2)

Screw-retained single crown in the molar region

ICA-Klebebasis_GH3 REF 31020
(ICA-Abutment GH3)

Screw-retained single crown in the molar region


With the scanbody Scan-Abutment Multi-Unit REF 31231:

MultiUnit_Locked_direkt
Konstr.m.Rotations
(direct construction
rotation-proof)

Construction directly on the Abutment Multi-Unit without Multi-Unit Cap titanium rotation-proof, e.g. a NEM (non-precious metal alloy) primary crown (milled in one piece, placed on the Multi-Units, screwed with the Screw Multi-Unit 3085M)

MultiUnit_NonLocked_direkt
Konstr.o.Rotations
(direct construction rotational)

Construction directly on the Abutment Multi-Unit without Multi-Unit-Cap titanium rotational, e.g. for bars (milled in one piece, directly placed on the Multi-Units, screwed with the Screw Multi-Unit 3085M)

MultiUnit_Cap Titan m. Rot.
REF 3111R
(titanium rotation-proof)

Constructions on Multi-Unit Cap titanium: e.g. screw-retained single crown, individual abutment (zircon), primary crown (zircon, NEM [non-precious metal alloy]), adhesives screwed with the Screw Multi-Unit 3085M

MultiUnit_Cap Titan o. Rot.
REF 3111
(titanium rotational)

Constructions on Multi-Unit Cap titanium: e.g. screw-retained bridge (zircon), individual bar (NEM [non-precious metal alloy]), adhesives screwed with the Screw Multi-Unit 3085M

2. Options for scanning Champion scanbodies

Intraoral Scanner/Laboratory Scanner

- 1) **Scanbody Impression Post Peek REF 3122 + Gingiva-Shuttle REF 3100**
Just like during impression taking, the Impression Post is placed intraorally on the Shuttle and scanned. In the case of a lab scanner, the Shuttle remains on the Champions (R)Evolution Laboratory Analog REF 31210; attach the Impression Post and scan.
- 2) **Scan-Abutment Champions (R)Evolution® REF 31230**
Remove the Shuttle from the implant either intraorally or on the model with the Laboratory Analog. Insert the Scan-Abutment, tighten with Screw 3085, and scan.
- 3) **Scan-Abutment Multi-Unit REF 31231**
 - a) First, screw the Multi-Unit Abutment (e.g., 0° GH1 REF 3108) intraorally or on the Champions (R)Evolution Laboratory Analog REF 31210 in the model. Then, place the Scan-Abutment Multi-Unit onto the Multi-Unit Abutment and tighten with the Screw Multi-Unit 3085M.
 - b) After taking the impression of the Multi-Units intraorally with the Impression Post-Set Multi-Unit REF 31223 and fabricating the model with the Multi-Unit Laboratory Analog REF 3114, tighten the Scan-Abutment Multi-Unit with the Screw Multi-Unit 3085M.

During the scanning process in the scanning software, follow the instructions on the screen for scanning a scanbody. Only place and scan the scanbody after being prompted by the software.

3. Instructions for selecting exocad

3a) Installation/Instructions exocad main library

It is recommended to have the installation performed by the exocad Reseller.


1. Download the file „Champions_Revolution_Exocad.zip“ from the Champions-Implants website at Download/CAD/CAM Library. Extract the file.
2. Copy the folder “Champions_Revolution” to the directory “DentalCADApp/library/implant/”.
3. If necessary, remove older versions of the Champions Revolution Library.
4. Restart exocad – from now on, the Champions implant system library is available.

3b) Installation exocad additional library Model Creator

1. Download the file “Zusatzbibliothek_exocad“ from the Champions-Implants website at Download/CAD/CAM Library. Extract the file.
2. Right-click with the mouse and select “Alle extrahieren“ (extract all) (confirm any warning dialog if prompted).
3. Copy the folders “Champions_Revolution_DIM“ and “Champions_Revolution_LabAnalogR“ to the directory “DentalCADApp/library/modelcreator/implants“.
4. If necessary, remove older versions of the Champions Revolution Library.
5. Restart exocad – from now on, the Champions implant system library is available.

3c) Selection options in exocad during construction


- a. Library selection: Champions – Revolution V2.14
- b. Scanbody selection: Abdruckpfosten REF 3122 (Impression Post Peek)
Scanabutment REF 31230
(Scan-Abutment Champions (R)Evolution®)
- c. Selection of construction options
From exocad Galway version onwards, the Abutment Adhesive 3104 in the hex can be rotated 6 times in the next step.


Selection 1st menu


Selection 2nd menu during the construction

Selection 3rd menu Impression Post

Selection 3rd menu Scan-Abutment

4. Installation/Instructions for selecting 3shape


It is recommended to have the installation of the Champions-library installed by the 3shape service partner.

For constructions based on Multi-Unit, extend the screw channel by 0.37 mm.


Selection options in 3shape in the order menu:

- a. Category selection: Champions Implants
- b. System selection: Champions – Revolution V2.14.1
Champions – Revolution V2.14.2
- c. Kit selection: Construction options
- d. Scanbody selection: Abdruckpfosten REF 3122 (Impression Post Peek)
Scanabutment REF 31230
(Scan-Abutment Champions (R)Evolution®)


Category selection


System selection

Kit selection

Scanbody selection

5. Instructions for Laboratory Analogs for 3D-printed models

REF 31210 Laboratory Analog Champions (R)Evolution

REF 31219 DIM-Analog

REF 31131 DIM-Analog Multi Unit


Creating the 3D-printed models using a Model Design software and compatible Champions-library.

Available in:
3shape Model Builder
exocad model creator
BISS Promadent
SHERAeasy-model


Only insert the Laboratory Analog REF 31210 with 3 grooves into the printed cavity at the implant position of the 3D-model. If necessary, reduce the horizontal distance for better stability. For REF 31219 DIM-Analog and REF 31131 DIM-Analog Multi-Unit, unscrew the DIM-Analog screw with the Prosthetic Driver REF 3086 and insert the Analog from the top of the 3D-printed model into the printed cavity at the implant position. Tighten the Analog screw from the bottom.

Selection of Laboratory Analog in exocad model creator

After (aligning the models) window name, the selection window appears.
Show model analog retention and make the model thicker.



Recommended model parameter settings:


Selection of Laboratory Analog in 3shape modelbuilder

The Laboratory Analog is determined when selecting the system.


REF 31219 DIM-Analog
REF 31131 DIM-Analog Multi-Unit

REF 31210 Laboratory Analog
Champions (R)Evolution

Recommended model parameter settings:
 Protect Laboratory Analog in the model with
 a chimney shape.

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Champions® is a registered trademark of
Champions-Implants GmbH
der Champions-Implants GmbH
de Champions-Implants GmbH