# PRODUKTINFORMATION

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## CONTENT

1. Generalities
2. Specifics
3. Immediate Restoration and Immediate Loading
4. Condensers and MIMI I-Technique
5. Drilling Sequences

### 1. Generalities

- 1) Please observe conservative basic rules and individual medical history of the patient.
- 2) Before an implant-supported and prosthetic restoration, among others periodontology and the overall situation should be considered and assessed.
- 3) Please follow the recommendations of the "Konsenskonferenz Implantologie" (consensus conference Implantology):  
  
In the upper jaw, a minimum of 6 implants/teeth is recommended for a removable denture and 8 implants/teeth for a fixed denture.  
  
In the lower jaw, a minimum of 4 implants/teeth is recommended for a removable denture and 6 implants/teeth for a fixed denture.
- 4) Pay attention to the sagittal jaw relation.
- 5) "Backward planning" is often advisable in many cases.
- 6) A tension-free fit of the superstructure ("passive fit") is essential for implant-supported dentures.
- 7) For the restoration of a single molar, a minimum of one  $\varnothing$  4.0 mm implant should be inserted.

## 2. Specifics

- 1) All Champions® implants can be inserted both using the minimally invasive method (without raising a mucoperiosteal flap, i.e. in a “flapless” way) and in a “classic” manner (with “full-flap”). They are particularly well suited for insertion following the secured MIMI® technique. An incision made with a scalpel or diamond turbine without raising a flap is also considered minimally invasive, as is a punch or direct penetration of the gingiva.
- 2) With Champions (R)Evolution® implants, the Shuttle (or the Gingiva-Clix to be clipped over it) should ideally terminate in a subgingival or equigingival position to completely exclude lateral micro-movements by the denture or the tongue during the first 6 to 8 weeks! If this requirement is not feasible, the Shuttle should be removed, and a surgical closure Hybrid Screw “Georgi” with Gingival Heights of 0.5 mm – 1.5 mm – 2.5 mm (instead of 3.5 mm with Shuttle) should be mounted hand-tightly.
- 3) Crestal relief in D1 & D2: Primary stability of 20 to 40 Ncm should be achieved in D1 and D2 bone exclusively through cancellous bone! Note that for example, with a ø 3.5 mm Champions (R)Evolution® or “New Art” implant, orange (crestal Drill ø 3.7 mm) and green (ø 4.0 mm) Drills should be used. All Champions® implants should ultimately have a primary stability of 20 to 40 Ncm, which corresponds roughly to hand-tight.
- 4) Pay attention to the settings of the Torque Wrench: Always start with the setting at 20 Ncm. Once the joint arm collapses, you have reached 20 Ncm. Now turn the middle mark from 20 to 40 Ncm with the handwheel. If it collapses again, you have then reached 40 Ncm, and so on.

### 3. Immediate Restoration and Immediate Loading

#### Immediate loading is usually possible when...

- ... a passive-fit denture is supported by a minimum of 4 implants/teeth.
- ... primary stability of 20 to 40 Ncm is achieved.
- ... with fewer than 4 implants/teeth (including single-tooth gaps), transitioning to the Secondary Osseointegration Stability phase (SOS) is recommended.  
The single crown should be integrated 8 weeks post-surgery!


#### Please note:

Champions® implants must only be placed and treated with the designated original Champions® instruments such as Drills, Condensers, Insertion Aids, and Prosthetic Drivers!


#### 4. Condensers and MIMI I-Technique

If the bone at the surgical site is soft (D3/D4), it is possible to transform soft bone into dense bone using Condensers in a minimally invasive way, which means to increase the D3/D4 bone densification, thus preparing it, for example, for immediate implantation. This process is known as "Osseous Metamorphosis." Additionally, Condensers assist implantologists in selecting the appropriate implant diameter because it's not the anatomy that determines the implant diameter, but rather the diameter of the corresponding Condenser that achieves primary stability. To do this, the implant bed is prepared using Condensers of increasing diameters.


- 1) If hand-tight stability and thus a primary stability of approximately 20 Ncm are achieved with the **ø 3.3 mm (long) Condenser**, a **3.5 mm diameter implant** can be inserted. If not, the cancellous bone is further compacted with the ø 3.8 mm Condenser.


- 2) If hand-tight stability and thus a primary stability of approximately 20 Ncm are achieved with the **ø 3.8 mm Condenser**, a **4.0 mm diameter implant** is used. If not, the cancellous bone is further prepared with the ø 4.3 mm Condenser.


- 3) If now hand-tight stability and thus a primary stability of approximately 20 Ncm are achieved with the **ø 4.3 mm Condenser**, a **4.5 mm diameter implant** is used. If not, the cancellous bone is further prepared with the ø 4.8 mm Condenser.

- 4) If now hand-tight stability and thus a primary stability of approximately 20 Ncm are achieved with the **ø 4.8 mm Condenser**, a **5.0 mm diameter implant** is inserted. If not, the cancellous bone is further prepared with the ø 5.3 mm Condenser.

- 5) If now hand-tight stability and thus a primary stability of approximately 20 Ncm are achieved with the **ø 5.3 mm Condenser**, a **5.5 mm diameter implant** is inserted.

## 5. Drilling Sequences


- 1) Whenever possible – both in the maxilla and mandible – start the pilot drilling 0.5 to 1 mm in a palatal or lingual position to the center of the alveolar ridge.
- 2) After anesthesia in the MIMI technique, measure the mucosal thickness with the yellow Drill (laser markings start at 2 mm, then 4 and 6 mm, and continue at 10, 12, 14 mm, etc.) until the instrument makes contact with bone.
- 3) Drilling should be done at a maximum of 250 rpm in the cortical bone and at 70 rpm in the cancellous bone.
- 4) To ensure primary stability in the cancellous bone of the D1 and D2 bone, it is generally recommended to provide “relief” of the crestal bone region by 1 to 2 mm vertically for titanium implants.
- 5) Typically, at least the full working length should be prepared for accommodating a Champions® implant, possibly with all Drills, preferably 1 to 2 mm longer if anatomical conditions allow and subcrestal placement is desired.
- 6) For D3 or D4 bone (e.g., maxillary posterior area), possibly prepare 2 mm less bone than the length of the implant so that an apical “grip” of 40 Ncm is achieved even with  $\varnothing$  3.0 mm or  $\varnothing$  3.5 mm Champions® implants. Also, during an Internal Direct Sinus Lift (IDS), drill only up to the “compact zone,” then use Condensers starting from  $\varnothing$  3.0 mm for sinus floor augmentation (lift) for 2 to 3 mm.

- 7) Cooling with sterile saline solution is not necessary. Especially with “immediate implants” (extraction and implantation in a single session), bone particles that are produced by the drilling should be “collected” (so you may be able to re-implant them later) or wiped off with a sterile gauze. By keeping the Drill clean, this increases cutting and compression capacity solely through drilling.
- 8) Never forget to perform the “BCC” (Bone Cavity Check) in all five dimensions (apical, buccal, lingual, mesial, distal) using the Champions® BCC Probe after the first and last drilling respectively.


Drilling and Condenser sequence for  $\varnothing$  3.5 mm  
Champions (R)Evolution® & Square "New Art" implants  
for all bone classes (D1 to D4):

- Pilot drilling with the yellow Drill to full working length.
- Preparation of the cavity with the black, white, and blue Drill respectively.
- Crestal relief with the brown, orange, and green Drill respectively.  
The green Drill should not be guided more deeply than 2 mm in a subcrestal position.
- Implant insertion

Drilling and Condenser sequence for  $\varnothing$  4.0 mm  
Champions (R)Evolution® & Square "New Art" implants  
for all bone classes (D1 to D4):

- Pilot drilling with the yellow Drill to full working length.
- Preparation of the cavity with the black, white, blue, and brown Drill respectively.
- Crestal relief with the orange, green, and red Drill respectively.  
The green and red Drills should not be guided more deeply than 2 mm in a subcrestal position.
- Implant insertion

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Champions® is a registered trademark of
Champions-Implants GmbH
der Champions-Implants GmbH
de Champions-Implants GmbH