

Ceramic implant now available

In releasing its new ceramic implants BioWin!, available as both one-piece and two-piece versions, Champions-Implants is breaking new ground. Without having a screw connecting the implant body and the abutment, the two-piece option is entirely metal-free. Since the glass fibre abutment is being glued, there is no space that is vulnerable to paropathogenic germs or bacteria, as it might happen in two-piece implant systems. The roughness of the surface is created using a one-of-a-kind patented process resulting in a faster osseointegration. Scientific studies (Prof. Becker, University Düsseldorf, among others) find osseointegration to be at 95.8 per cent. The implants are available in five different lengths (6, 8, 9, 11, 13 mm) and diameters (4.1, 4.5, 5 mm). Moreover, there is a free and user-friendly software available, with which individual implants having individual emergence profiles can be developed. BioWin! implants can be inserted using either a minimally invasive approach or a classic full-flap one.



Champions-Implants GmbH
Champions Platz 1
55237 Flonheim, Germany
www.championsimplants.com